

Specially made:

Vena® Vario Levant

The Vena® Vario Levant phlebotomy chair combines comfort and seamless functionality. The easy operation of Vena® Vario Levant makes it possible to act quickly in emergency situations. The patient can be put in the Trendelenburg position in one quick move.

Vena® Vario Levant makes it possible to sit comfortably and relax during blood sampling. The nurse can provide back and leg support simultaneously which ensures good support to the patient's whole body. When reclined, the nurse or medical expert can also use Vena® Vario Levant as treatment table. The height is adjustable by foot pump.

Vena® Vario Levant is equipped with two height adjustable armrests for the left and right arm. By opting for the multi-adjustable armrests it is even easier to put the patient's arm in a correct and relaxed position and the nurse can draw blood from the patient's arm in an ergonomically responsible manner.

Specifications:

- Comfortable high backrest for relaxing and sitting comfortably
- Choose PUxx upholstery for cleaning with alcohol (up to 70%) and for easy cleaning of (blood) stains
- Trendelenburg position in 2 seconds!
- Get-up brackets to facilitate easy sitting down and standing up
- Simple design prevents the accumulation of dirt and makes cleaning easy
- Choose from different functional and ergonomic armrests
- Hospital grade swivel castors for manoeuverability and transport
- Equipped with hydraulic foot pump for 20 cm height adjustment, even in Trendelenburg position
- Extra ease of use by choosing a central brake on both front castors
- Optional foot operated back- and leg rest adjustment









Standard model:

- Armrests: 2 x 'Vena® fxed' Pur
- Height adjustment: 47 67 cm
- Sizes.: reclined position167,5 cm, base 58 (W) cm x 70 (L) cm
- Hospital grade castors (with brake)
- Robust frame (white RAL 9016)
- Trendelenburg position quick release

Options/accessories:

- Armrests: 'Vena® multi adjustable' in PU or upholstered in Stamskin or PUxx
- Get-up brackets
- Central brake on both front castors
- Foot operated back- and leg rest adjustment

Upholstery

Stamskin upholstery:

Stamskin® 100% waterproof; Urine-, blood-, sweat-, fungal proof. Clean with water and soap, rinse of with water.



PUxx upholstery:

PUxx® 100% waterproof, PVC and Phthalate free; resistant against commonly used cleaning products and disinfectants (alcohol based disinfectants, active chlorine, quaternary ammonium compounds).

